# Longsee-Antigentest

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Stand 14.12.2021

# DECLARATION OF CONFORMITY

**Manufacturer** Guangdong Longsee Biomedical Co.,Ltd.  
**Address** 5/F Building A, No.83, Ruihe Road, Huangpu District, 510000, Guangzhou, China

**European Representative** MedPath GmbH  
**Address** Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

**Product Information** 2019-nCoV Ag Rapid Detection Kit (Immuno-Chromatography)  
**Model code** LS-C-T-009  
**Classification** Other IVD Device  
**Registration Number in German DIMDI Database** DE/CA61/1M50/294  
**Conformity Assessment Route: Annex III**

**General Applicable Directives:** In vitro diagnostic medical devices Directive: 98/79/EC

<b>Standards Applied</b>	EN13612:2002/AC: 2002	EN ISO13485:2016
	EN ISO 23640:2015	EN ISO14971:2012
	EN 13641:2002	EN ISO18113-1:2011
	EN 15223-1:2016	EN ISO18113-2:2011

We, the manufacturer, hereby declare under our sole responsibility that the above mentioned products meet the provisions of the following EC Council Directives and Standards. The products meet prospective uses and all supporting documentations are retained under the premises of the manufacturer.

Place, date of issued: Guangzhou, P. R. China, May 10, 2021

Signature of Vice President:

  
(Marco Zhang)

## 2019-nCoV Ag Rapid Detection Kit(Immuno-Chromatography)

Instructions for Use

LS-C-T-009

### 【Intended Use】

The kit is an immunochromatographic test for rapid, qualitative detection of the SARS-CoV-2 nucleocapsid antigen extracted from the nasopharyngeal swab or oropharyngeal swab specimens from persons suspected of COVID-19. The kit is intended to support the rapid diagnosis of SARS-CoV-2 infections.

**For in vitro diagnostic use only. For professional use only.**

### 【Summary】

The novel coronavirus belongs to the  $\beta$  genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

### 【Test Principle】

This kit applies immunochromatography technology to detect the presence or absence of 2019-nCoV nucleocapsid proteins in swab specimens from patients with signs and symptoms of infection suspected of having 2019-nCoV by the double antibody sandwich method. While the concentration of the 2019-nCoV antigens in samples is higher than or equal to the minimum detection limit, these antigens react separately with corresponding antibodies to form complexes, and the 2019-nCoV antibodies are coated in the detection area (T). These antigens are captured and a red line of reaction is formed. The result is rated as positive. Otherwise the result formed in T without a red line is assessed as negative. Under normal test conditions, the quality control area (C) should be colored to indicate that the test is valid.

### 【Component】

	Component	25 T/kit	50 T/kit	100 T/kit	200 T/kit
1	Test Cartridge	25 pcs	50 pcs	100 pcs	200 pcs
2	One-off Swab	25 pcs	50 pcs	100 pcs	200 pcs
3	Extraction Tube	25 pcs	50 pcs	100 pcs	200 pcs

### 【Storage and Stability】

- 1) Store at 4-35 °C up to the expiration date printed on the package.
- 2) Do not freeze the kit or its components.

### 【Specimen Collection and Preparation】

The test can be performed with oropharyngeal swab or nasopharyngeal swab specimen.

1. Oropharyngeal swab specimen collection: Insert swab into the posterior pharynx and tonsillar areas. Rub swab over both tonsillar pillars and posterior oropharynx and avoid touching the tongue, teeth, and gums. Place swab, tip first, into the extraction tube provided.
2. Nasopharyngeal swab specimen collection: Tilt patient's head back 70 degrees. Gently and slowly insert a minitip swab with a flexible shaft (wire or plastic) through the nostril parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx. Gently rub and roll the swab. Leave swab in place for several seconds to absorb secretions. Slowly remove swab while rotating it. Specimens can be collected from both sides using the same swab, but it is not necessary to collect specimens from both sides if the minitip is saturated with fluid from the first collection. If a deviated septum or blockage create difficulty in obtaining the specimen from one nostril, use the same swab to obtain the specimen from the other nostril. Place swab, tip first, into the extraction tube provided.
3. It is recommended that specimens be treated with the sample extraction buffer provided with the kit as soon as possible after collection. If immediate processing is not possible, the specimen can be stored in a dry, sterilized and tightly sealed plastic tube at 2°C ~8°C for up to 8 hours.

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#### Test Preparation


1. Allow all kit components to equilibrate to room temperature (15-30°C) prior to testing for 30 minutes, if previously stored in a cool place.

### Extraction

1. Open the drip cap of the extraction tube and immerse the sampled swab tip into the extraction tube and rotate the swab tip 10 times in the buffer liquid while applying pressure with your fingers. Then let it sit at room temperature for 5 min.
2. Remove the swab while squeezing the sides of the tube to extract the liquid with the sample from the swab.
3. Screw the drip cap tightly onto the extraction tube.
4. Screw off the tip of the drip cap of the extraction tube.

### Reaction with Test Cartridge

5. Remove a test cartridge from the sealed pouch by tearing at the notch and place it on a level surface. Drip 3-4 drops (about 100  $\mu$ L) of extracted specimens vertically into the specimen well (S) on the test cartridge by squeezing the tube. The formation of air bubbles in the specimen well (S) must be avoided. Do not handle or move the test cartridge until the test is complete and ready for reading.
6. Start timer. Read result within 15-20 minutes of adding the sample. The test result is invalid after 20 minutes.

### 【Interpretation of the Result】

To read the test results, all you have to do is look at the results window.

1. Positive


If two color bands respectively appear at the position of the quality control line (C line) and the detection line (T line) in the observation window, the test result is positive. The test result indicates that the sample contains 2019-nCoV antigens.

2. Negative

If one color band appear at the position of the quality control line (C line) and no one color band at the detection line (T line) in the observation window, the test result is positive negative. The test result indicates that the 2019-nCoV antigens in the sample are negative or the concentration is below the detection limit of the set.

3. Invalid

If no color band appear at the position of the quality control line (Line C) or only one color band at the detection line (T line) in the observation window, the test is invalid. The sample should be taken again and the test repeated.

### 【Limitations】

1. The contents of this kit are for professional use and qualitative detection of 2019-nCoV antigen from swab specimen. Other specimen types may lead to incorrect results and must not be used.
2. Failure to follow the instructions for test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
3. Due to the limitations of the methodology, experimenters should pay more attention to negative results. A negative test result may occur if the specimen was collected, extracted or transported improperly. A negative test result does not eliminate the possibility of 2019-nCoV infection and should be confirmed by viral culture or a molecular assay.
4. Positive test results do not rule out co-infections with other pathogens. Positive results may occur in cases of infection with SARS-CoV.
5. Test results must be evaluated in conjunction with other clinical data available to the physician.
6. Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.

### 【Performance Characteristics】

1. LIMIT OF DETECTION

The limit of detection has been evaluated at  $6.0 \times 10^2$  TCID<sub>50</sub>/mL.

2. CROSS-REACTIVITY

There was no cross-reaction and interference with the potential cross-reacting microorganisms listed below:

	Potential Cross-Reactant	Test Concentration
Virus	Adenovirus	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Human coronavirus 229E	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Human coronavirus OC43	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Human coronavirus NL63	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Human coronavirus HKU1	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	MERS-coronavirus	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	SARS-coronavirus	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Human Metapneumovirus (hMPV)	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Parainfluenza virus 1	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Parainfluenza virus 2	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Parainfluenza virus 3	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Parainfluenza virus 4	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Influenza A	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Influenza B	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Enterovirus	1.0×10 <sup>7</sup> TCID <sub>50</sub> /mL
	Respiratory syncytial virus	1.0×10 <sup>7</sup> PFU/mL
Rhinovirus	1.0×10 <sup>8</sup> IFU/mL	
Bacteria	Bordetella pertussis	1.0×10 <sup>6</sup> cells/mL
	Chlamydia pneumoniae	1.0×10 <sup>6</sup> IFU/mL
	Haemophilus influenzae	1.0×10 <sup>6</sup> cells/mL
	Legionella pneumophila	1.0×10 <sup>6</sup> cells/mL
	Mycoplasma pneumoniae	1.0×10 <sup>6</sup> U/mL
	Streptococcus pyogenes	1.0×10 <sup>6</sup> cells/mL
	Streptococcus pneumoniae	1.0×10 <sup>6</sup> cells/mL
	Mycobacterium tuberculosis	1.0×10 <sup>6</sup> cells/mL
	Staphylococcus aureus	1.0×10 <sup>6</sup> org/mL
	Staphylococcus epidermidis	1.0×10 <sup>6</sup> org/mL
Yeast	Candida albicans	1.0×10 <sup>6</sup> cells/mL

### 3. INTERFERING SUBSTANCES

There was no interference for potential interfering substances listed below:

Substance	Concentration
Whole Blood	4%
Mucin	0.5%
Chloraseptic (Menthol/Benzocaine)	1.5 mg/mL
Naso GEL (NeilMed)	5% v/v
CVS Nasal Drops (Phenylephrine)	15% v/v
Afin (Oxymetazoline)	15% v/v
CVS Nasal Spray (Cromolyn)	15% v/v
Zicam	5% v/v
Homeopathic (Alkalol)	1:10 dilution
Sore Throat Phenol Spray	15% v/v
Tobramycin	4 µg/mL
Mupirocin	10 mg/mL
Fluticasone Propionate	5% v/v
Tamiflu (Osetamivir Phosphate)	5 mg/mL

### 4. HOOK EFFECT

There is no hook effect at 1.5×10<sup>7</sup>TCID<sub>50</sub>/mL of SARS-CoV-2 which was isolated from a COVID-19 confirmed patient in China.

### 5. CLINICAL EVALUATION

Clinical evaluation was performed to compare the results obtained by The SARS-CoV-2 Antigen Rapid Test Kit and RT-PCR. The results were summarized below:

#### 5.1 Nasopharyngeal swabs

Nasopharyngeal Swabs	Reference Extracted SARS-CoV-2 RT-PCR assay				95% CI			
	POS	NEG	Total	PPA	LCI	UCI		
2019-nCoV Ag Rapid Detection Kit (Immuno-Chromatography)	POS	340	1	341	NPA	99.72%	98.23%	99.99%
	NEG	16*	362	378	PPV	99.71%	98.12%	99.98%
	Total	356	363	719	NPV	95.77%	93.08%	97.48%
				Prevalence	49.51%	45.80%	53.23%	
				OPA	97.64%	96.25%	98.52%	
				Kappa	95.27%	93.43%	97.11%	

\*3 of the discrepant samples had high Ct values (>30) when tested by the comparative method

#### 5.2 Oropharyngeal swabs

Oropharyngeal Swabs	Reference Extracted SARS-CoV-2 RT-PCR assay				95% CI			
	POS	NEG	Total	PPA	LCI	UCI		
2019-nCoV Ag Rapid Detection Kit (Immuno-Chromatography)	POS	339	1	340	NDA	99.72%	98.23%	99.99%
	NEG	17*	362	379	PPV	99.71%	98.11%	99.98%
	Total	356	363	719	NPV	95.51%	92.78%	97.28%
				Prevalence	49.51%	45.80%	53.23%	
				OPA	97.50%	96.08%	98.41%	
				Kappa	94.99%	93.10%	96.88%	

\*4 of the discrepant samples had high Ct values (>30) when tested by the comparative method

#### Marknote:

PPA - Positive Percent Agreement (Sensitivity)  
NPA - Negative Percent Agreement (Specificity)  
PPV - Positive Predictive Value

NPV - Negative Predictive Value  
OPA - Overall Percent Agreement  
CI - Confidence Interval

LCI - Lower Confidence Interval  
UCI - Upper Confidence Interval

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Symbols	Title of symbol	Symbols	Title of symbol	Symbols	Title of symbol
	In vitro diagnostic device		Do not re-use		manufacturer
	Avoid direct sunlight		Use-by date		Batch code
	Consult instructions for use		Do not use if package is damaged		Temperature limit
	CE Mark		Authorized representative in the European Community		

### 【Warnings and Precautions】

1. Take the appropriate precautions for infected material if necessary.
2. There is a desiccant in the package and it is not allowed to be taken orally.
3. Please read the instructions carefully before use and strictly follow the instructions. The reagents cannot be used when expired or the packaging bag is damaged or the seal fails. There should be no coloured lines on the test card before use.
4. When opening the foil pouch, pay attention that the test cartridge should not fall out. After the foil pouch is opened, test cartridge should be used within 30 minutes .
5. The test kit should be sealed and stored at 4-35 °C, keeping away from moisture and sunlight, the test kit stored at low temperature should be restored to room temperature before use.
6. For samples containing or suspected of containing the source of infection, there should be appropriate bio-safety guarantee safe operation procedures. The following are relevant precautions:
  - 1) The gloves or reagents used in processing samples should be disinfected,
  - 2) Use disinfectant to disinfect the spilled samples or reagents.
7. It is a disposable in vitro diagnostic product. For experimental wastes such as test cards, gloves, pipette tips, unused samples or reagents, etc., which have potential biological hazards, should be in accordance with biological safety regulations, environmental protection regulations or medical waste regulations for disinfection and disposal.
8. Avoid skin or eyes touching the sample extraction buffer. If accidentally touch sample extraction buffer, immediately use plenty water to wash the skin or eyes. Forbid drinking the sample extraction buffer. If carelessly intake, please gargle thoroughly. Go to hospital if you feel not well.

**Guangdong Longsee Biomedical Co.,Ltd.**  
5/F Building A, No.83, Ruihe Road, Huangpu District, 510000, Guangzhou, China  
Contact number: + 86-20-66234660  
Website: <http://www.longseemed.com/>

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**Guangdong Longsee Biomedical Co., Ltd.**